

---

**flexlib**

***Release 1.1.1***

**Jonas Walkling (MinkTec)**

**Aug 01, 2025**

## CONTENTS:

<b>1</b>	<b>Coordinate Systems and Orientation</b>	<b>2</b>
<b>2</b>	<b>Measurement</b>	<b>5</b>
<b>3</b>	<b>SensorAngles</b>	<b>9</b>
<b>4</b>	<b>Coordinates</b>	<b>11</b>
<b>5</b>	<b>CaseOrientation</b>	<b>12</b>
<b>6</b>	<b>FullFlexTailReconstructionPy</b>	<b>13</b>
<b>7</b>	<b>Common Metrics and Analysis</b>	<b>14</b>
<b>8</b>	<b>Distance Parameter Usage Examples</b>	<b>15</b>
<b>9</b>	<b>flexlib</b>	<b>17</b>
9.1	flexlib package . . . . .	17
<b>10</b>	<b>Demo: Plotting</b>	<b>33</b>
<b>11</b>	<b>Dishwasher Demo</b>	<b>39</b>
11.1	Schmitt Trigger for Bending Detection . . . . .	41
11.2	Summary . . . . .	43

**flexlib** is a Python library for processing and analyzing sensor data from FlexTail wearable sensors. The library provides tools for loading, processing, and evaluating biomechanical measurements.

**Core Functionality:**

- **Data Processing:** Load and process sensor measurements from various file formats (RSF v1/v2, CSV)
- **Biomechanical Analysis:** Calculate spine angles, posture metrics, and movement parameters
- **Signal Processing:** Apply filters for motion and bending detection.
- **Data Evaluation:** Extract metrics from measurement data for analysis

**Key Components:**

- Measurement objects representing individual sensor readings with timestamp, flex sensor data, accelerometer, and gyroscope values
- AnnotatedRecording for time-series data with event annotations
- Evaluation metrics for extracting biomechanical features (lumbar angles, twist, lateral/sagittal flexion)
- Signal processing utilities including Schmitt triggers for robust motion detection
- File I/O support for reading and writing sensor data in multiple formats

## COORDINATE SYSTEMS AND ORIENTATION

For anatomical reference, the following orientation planes are used:




Figure 1.1: **Anatomical Planes:** Reference planes for understanding sensor orientation relative to human anatomy.

All references to sagittal or lateral are in reference to these planes.

The sensor comes in two variants: 36 cm and 45 cm. The amount of sensor pairs is the same which means that the distance between the measuring points is 20mm for the 36 cm variant and 25mm for the 45 cm variant. To handle different sensor variants, specify the appropriate `dist` parameter when creating a `Measurement` object:

- For 36cm sensors: `flexlib.Measurement(data, time, dist=20.0)` (default)
- For 45cm sensors: `flexlib.Measurement(data, time, dist=25.0)`

The distance parameter can also be overridden for individual calculations using `coordinates_with_dist()` or `full_reconstruction_with_dist()` methods, allowing you to compare results with different sensor spacings without recreating the measurement object.




Figure 1.2: **Sensor Overview:** A drawing of the sensor strip and the case including directions of the accelerometer and gyroscope. The `acc`, `gyro`, `left` and `right` are lists of values that are used in the measurement class.

The main sensor data fields are:

- **left:** List of integer values from the left side of the sensor strip. Positive values indicate bending forward, negative values indicate bending backward.
- **right:** List of integer values from the right side of the sensor strip. Positive values indicate bending forward, negative values indicate bending backward.
- **gyro:** Gyroscope readings. This is a list of three integer values representing angular velocity (rotation rate) around the X, Y, and Z axes. Used to analyze rotational motion.
- **acc:** Accelerometer readings. This is a list of three integer values representing acceleration along the X, Y, and Z axes. A value of 2000 corresponds to 1g (gravitational acceleration). Used to determine orientation and movement.




Figure 1.3: **Coordinate System:** All coordinates that are calculated by the sensor follow the coordinate system shown here and are in mm. The X axis is the right/left axis in the dorsal plane, the Y axis is the forward/backward axis in the midsagittal plane, and the Z axis is the vertical axis.




Figure 1.4: **Case Orientation:** We measure the case orientation by two angles: pitch and roll. Those angles are calculated from the accelerometer data and represent the tilt of the sensor case in space. Pitch is positive to the front, roll is positive to the right.

## MEASUREMENT

**class** flexlib.**Measurement**(*data, time=None, dist=20.0*)

The Measurement class is the core representation of a FlexTail measurement. It encapsulates the raw sensor data (sensor strip, accelerometer and gyroscope), the number of sensors, and the timestamp of the measurement.

### Parameters

- **data** (*list[int]*) – Raw sensor, accelerometer, and gyroscope data
- **time** (*int or None*) – Timestamp in ms since epoch
- **dist** (*float*) – Distance between sensor elements in mm (default: 20.0 for 36cm sensors, use 25.0 for 45cm sensors)

### values

Raw sensor data values.

#### Type

list[int]

### timestamp

Timestamp in milliseconds since epoch.

#### Type

int or None

### dist

Distance between sensor elements in mm.

#### Type

float

### datetime

Python datetime object from timestamp.

#### Type

datetime.datetime or None

### iso\_time

ISO formatted time string.

#### Type

str

### left

The measured values from the left sensor side. Positive values indicate bending forward, negative values indicate bending backward.

#### Type

list[int]

**right**

The measured values from the right sensor side. Positive values indicate bending forward, negative values indicate bending backward.

**Type**

list[int]

**acc**

Accelerometer readings. Value of 2000 corresponds to 1g (gravitational acceleration).

**Type**

list[int]

**gyro**

Gyroscope readings (angular velocity around each axis).

**Type**

list[int]

**voltage**

Voltage reading.

**Type**

int

**angles**

Angles between adjacent measurement points along the sensor strip.

**Type**

*SensorAngles*

**coordinates**

3D coordinates of each measurement point along the sensor strip, calculated from the bending angles.

**Type**

*Coordinates*

**case\_orientation**

Orientation of the sensor case, derived from accelerometer data.

**Type**

*CaseOrientation*

**number\_of\_sensors**

Number of sensors in the strip.

**Type**

int

**lateral\_flexion**

Lateral flexion angle. Measures lateral (side-to-side) flexion. Useful for analyzing sideways bending.

**Type**

float

**sagittal\_flexion**

Sagittal flexion angle. Measures sagittal (forward/backward) flexion. This is the angle between the case and the tip of the sensor. If person would lean to the right the angle would be positive. Useful for analyzing bending in the sagittal plane.

**Type**

float



**coordinates\_with\_dist(*dist*)**

3D coordinates of each measurement point along the sensor strip with custom distance parameter.

**Parameters**

**dist** (*float* or *None*) – Distance between sensor elements in mm. If *None*, uses the measurement's dist value.

**Returns**

3D positions in mm

**Return type**

*Coordinates*

**full\_reconstruction()**

Get full 3D reconstruction of the sensor strip.

**Returns**

Full reconstruction data

**Return type**

*FullFlexTailReconstructionPy*

**full\_reconstruction\_with\_dist(*dist*)**

Get full 3D reconstruction of the sensor strip with custom distance parameter.

**Parameters**

**dist** (*float* or *None*) – Distance between sensor elements in mm. If *None*, uses the measurement's dist value.

**Returns**

Full reconstruction data

**Return type**

*FullFlexTailReconstructionPy*

**to\_csv\_row()**

Convert measurement to CSV row format.

**Returns**

CSV formatted string

**Return type**

str

**movement()**

Calculate movement metric from gyroscope data.

**Returns**

Movement value

**Return type**

float

**calc\_sagittal\_approx()**

Calculates the sagittal approximation based on the angles. The approximation assumes that the sensor does not deform along the x axis.

**Returns**

Sagittal approximation

**Return type**

float

**calc\_lateral\_approx()**

Calculate lateral approximation.

**Returns**

Lateral approximation

**Return type**

float

## SENSORANGLES

```
class flexlib.SensorAngles(left, right)
    Bending and rotational angles between measurement points.
```

- Measuring point
- Sensor strip (viewed from the right)
- - - Connection between adjacent points
- ⤿ Angle that is measured




Figure 3.1: **Sensor Angles:** Visualization of alpha (bending) angles between adjacent points along the sensor strip. The black squares represent measurement points on the sensor. The numbers are the readouts produced by each specific measuring point. Each point measures the angle between the straight line connecting it (n) to the next point (n+1) and the connection between the next (n+1) and the point after that (n+2). The angles depicted here are the alpha values. The beta values represent the rotation or twist between adjacent measurement points.

### Parameters

- **left** (*list[int]*) – Left sensor readings
- **right** (*list[int]*) – Right sensor readings

### alpha

Bending angles (radians) between adjacent points (sagittal plane).

#### Type

*list[float]*

### beta

Rotational angles (radians) between adjacent points (twisting).

#### Type

*list[float]*

**bend**

The overall bending (flexion/extension) of the lumbar region, calculated as the sum of the bottom 9 alpha values. Useful for posture and spinal movement analysis.

**Type**

float

**twist**

The total rotational movement (twisting) of the sensor strip, calculated as the sum of all beta angles. Useful for detecting torsional motions of the back.

**Type**

float

**static calc\_angles**(*left*, *right*)

Calculate angles from left and right sensor data using default distance (20.0mm).

**Parameters**

- **left** (*list[int]*) – Left sensor readings
- **right** (*list[int]*) – Right sensor readings

**Returns**

Tuple of (alpha, beta) angles

**Return type**

tuple[list[float], list[float]]

**static calc\_angles\_with\_dist**(*left*, *right*, *dist*)

Calculate angles from left and right sensor data with custom distance parameter.

**Parameters**

- **left** (*list[int]*) – Left sensor readings
- **right** (*list[int]*) – Right sensor readings
- **dist** (*float*) – Distance between sensor elements in mm

**Returns**

Tuple of (alpha, beta) angles

**Return type**

tuple[list[float], list[float]]

**static new\_with\_dist**(*left*, *right*, *dist*)

Create SensorAngles with custom distance between sensor elements.

**Parameters**

- **left** (*list[int]*) – Left sensor readings
- **right** (*list[int]*) – Right sensor readings
- **dist** (*float*) – Distance between sensor elements in mm

**Returns**

SensorAngles object with calculated angles

**Return type**

*SensorAngles*

## COORDINATES

**class** flexlib.Coordinates(*x*, *y*, *z*)

3D positions of each measurement point along the sensor strip. The coordinate system is shown in [Figure 1.3](#).

**Parameters**

- **x** (*list[float]*) – X-coordinates (right/left in dorsal plane, mm)
- **y** (*list[float]*) – Y-coordinates (forward/backward in midsagittal plane, mm)
- **z** (*list[float]*) – Z-coordinates (vertical, mm)

**x**

X-coordinates (right/left in dorsal plane, mm).

**Type**

*list[float]*

**y**

Y-coordinates (forward/backward in midsagittal plane, mm).

**Type**

*list[float]*

**z**

Z-coordinates (vertical, mm).

**Type**

*list[float]*

**rotated**(*pitch*, *roll*)

Returns a new Coordinates with all points rotated by the given pitch (x-axis) and roll (y-axis) angles (in radians).

**Parameters**

- **pitch** (*float*) – Rotation around x-axis (radians)
- **roll** (*float*) – Rotation around y-axis (radians)

**Returns**

Rotated coordinates

**Return type**

*Coordinates*

## CASEORIENTATION

**class** flexlib.**CaseOrientation**(*pitch*, *roll*)

Orientation of the sensor case in space.

**Parameters**

- **pitch** (*float*) – Angle in sagittal plane (forward/backward tilt, radians)
- **roll** (*float*) – Angle in dorsal plane (side-to-side tilt, radians)

**pitch**

Angle in sagittal plane (forward/backward tilt, radians).

**Type**

float

**roll**

Angle in dorsal plane (side-to-side tilt, radians).

**Type**

float

**static from\_acc**(*acc*)

Calculate case orientation from accelerometer data.

**Parameters**

**acc** (*list[int]*) – Accelerometer readings

**Returns**

Case orientation

**Return type**

*CaseOrientation*

## FULLFLEXTAILRECONSTRUCTIONPY

### **class** flexlib.FullFlexTailReconstructionPy

3D reconstruction of the sensor strip, including all calculated coordinates and angles. The `left` and `right` coordinates represent the left / right outline of the sensor strip, when viwed from the top. The `center` coordinates represent the center of the sensor strip, and are equivalent to the `coords` attribute.

#### **alpha**

Bending angles along the midsagittal plane (forward/backward flexion).

##### **Type**

`list[float]`

#### **beta**

Rotational angles (twisting) along the sensor strip.

##### **Type**

`list[float]`

#### **coords**

3D coordinates of the sensor strip in space.

##### **Type**

*Coordinates*

#### **case\_position**

Orientation of the sensor case (pitch, roll).

##### **Type**

*CaseOrientation*

#### **left**

3D coordinates for left points of the sensor tip.

##### **Type**

*Coordinates*

#### **center**

3D coordinates for center points of the sensor tip.

##### **Type**

*Coordinates*

#### **right**

3D coordinates for right points of the sensor tip.

##### **Type**

*Coordinates*

## COMMON METRICS AND ANALYSIS

The FlexTail sensor system provides several key metrics for biomechanical analysis:

### Primary Measurements:

- **LUMBAR\_ANGLE (bend):** Measures the overall bending (flexion/extension) of the lumbar region. Calculated by summing the bottom 9 alpha values. Useful for posture and spinal movement analysis.
- **TWIST:** Measures rotational movement (twisting) of the sensor strip. Calculated by summing all beta angles. Useful for detecting torsional motions of the back.
- **LATERAL:** Measures lateral (side-to-side) flexion. Useful for analyzing sideways bending.
- **SAGITTAL:** Measures sagittal (forward/backward) flexion. Useful for analyzing bending in the sagittal plane.

### Approximation Methods:

- **LATERAL\_APPROX:** Provides an approximation of lateral flexion using a simplified calculation. Useful for quick or less precise assessments.
- **SAGITTAL\_APPROX:** Provides an approximation of sagittal flexion using a simplified calculation. Useful for quick or less precise assessments.

### Specialized Metrics:

- **THORACIC\_ANGLE:** Focuses on the bending angle in the thoracic (upper back) region, typically using the second half of the sensor's angles. Useful for upper back posture analysis.
- **ACCELERATION:** Extracts the acceleration along a specific axis (usually the main axis of the sensor). Useful for detecting movement or orientation changes.
- **GYRO:** Measures angular velocity (rotation rate) from the gyroscope. Useful for analyzing dynamic movements and rotations.



## DISTANCE PARAMETER USAGE EXAMPLES

The FlexTail library supports different sensor variants with varying element spacing. Here are practical examples of using the distance parameter functionality:

### Basic Usage with Different Sensor Types:

```
import flexlib

# Load measurement data
data = [100, 150, 200, ...] # Your sensor data

# For 36cm sensors (20mm spacing) - default
measurement_36cm = flexlib.Measurement(data, dist=20.0)

# For 45cm sensors (25mm spacing)
measurement_45cm = flexlib.Measurement(data, dist=25.0)
```

### Method-Level Distance Overrides:

```
# Create measurement with default 20mm spacing
measurement = flexlib.Measurement(data)

# Get coordinates with different spacing values
coords_default = measurement.coordinates # Uses 20.0mm
coords_15mm = measurement.coordinates_with_dist(15.0) # Override to 15mm
coords_25mm = measurement.coordinates_with_dist(25.0) # Override to 25mm

# Get full reconstruction with custom spacing
recon_default = measurement.full_reconstruction # Uses 20.0mm
recon_30mm = measurement.full_reconstruction_with_dist(30.0) # Override to 30mm
```

### Custom Angle Calculations:

```
# Calculate angles with specific sensor spacing
left_data = [100, 110, 120, ...]
right_data = [95, 105, 115, ...]

# Default 20mm spacing
angles_default = flexlib.SensorAngles(left_data, right_data)

# Custom 25mm spacing
angles_25mm = flexlib.SensorAngles.new_with_dist(left_data, right_data, 25.0)

# Raw angle calculation (returns tuple)
alpha, beta = flexlib.SensorAngles.calc_angles_with_dist(left_data, right_data, 30.0)
```

### Comparing Different Sensor Configurations:

```
# Compare coordinate scaling effects
measurement = flexlib.Measurement(data)

distances = [15.0, 20.0, 25.0, 30.0]
for dist in distances:
    coords = measurement.coordinates_with_dist(dist)
    last_point = coords.x[-1], coords.y[-1], coords.z[-1]
    print(f"Distance {dist}mm - Last point: {last_point}")
```

This approach ensures accurate physics-based calculations rather than post-processing scaling, providing more reliable biomechanical analysis results.

## 9.1 flexlib package

### class AbstractSchmittTrigger

Bases: Generic[T]

Base class for Schmitt triggers with hysteresis-based state detection. The Schmitt trigger classes implementations for robust signal state detection based on hysteresis. Schmitt triggers are fundamental components in signal processing and digital electronics, used to convert noisy or analog input signals into clean, discrete digital outputs. They achieve this by introducing hysteresis, which means the trigger has two different threshold levels: one for transitioning from LOW to HIGH (the upper threshold), and another for transitioning from HIGH to LOW (the lower threshold). This prevents rapid toggling or “chattering” of the output when the input signal is noisy or hovers near the threshold. This module provides several flexible and extensible Schmitt trigger implementations:

#### **\_state**

The current state of the trigger (HIGH or LOW).

#### **Type**

*SchmittState*

#### **add(x)**

Add new input value and return new state if changed.

#### **Parameters**

**x** (*T*) – Input value to evaluate.

#### **Returns**

New state if changed, None if no state change.

#### **Return type**

*SchmittState* or None

#### **property state: SchmittState**

Current trigger state.

#### **Returns**

The current state (HIGH or LOW).

#### **Return type**

*SchmittState*

### class AnnotatedRecording(measurements, annotations=[])

Bases: object

A collection of sensor measurements with time-based annotations.

Contains a time-series of sensor measurements and associated annotations that mark specific events or states during the recording period.

**property end\_time**

Get the timestamp of the last measurement.

**static from\_json\_entry**(*json*, *dist*=20.0)

Create an AnnotatedRecording from a JSON representation.

**Parameters**

- **json** – Dictionary with ‘annotations’ and ‘content’ keys
- **dist** (*float*) – Distance between sensor elements in mm (default: 20.0)

**Returns**

AnnotatedRecording instance

**Return type**

*AnnotatedRecording*

**split\_by\_annotations**()

Split the recording into segments based on annotations.

**Returns**

List of tuples containing (annotation\_label, measurements\_in\_segment)

**Return type**

*List*[*Tuple*[str, *List*[*Measurement*]]]

**property start\_time**

Get the timestamp of the first measurement.

**to\_json\_entry**()

Convert the annotated recording to a JSON-serializable dictionary.

**Returns**

Dictionary with ‘annotations’ and ‘content’ keys

**Return type**

dict

**class CSVReader**

Bases: object

Parse sensor measurement data from CSV files or strings.

**static parse**(*path*, *dist*=20.0)

Parse CSV content from a file.

**Parameters**

- **path** (*str*) – Path to the CSV file to parse.
- **dist** (*float*) – Distance between sensor elements in mm (default: 20.0)

**Returns**

Parsed measurements in an annotated recording.

**Return type**

*AnnotatedRecording*

**static parse\_string**(*content*, *dist*=20.0)

Parse CSV content from a string.

**Parameters**

- **content** (*str*) – The CSV content as a string.
- **dist** (*float*) – Distance between sensor elements in mm (default: 20.0)

**Returns**

Parsed measurements in an annotated recording.

**Return type**

*AnnotatedRecording*

**class CSVWriter**

Bases: object

A writer for exporting measurement data to CSV format.

Writes sensor values followed by timestamp in comma-separated format.

**static write**(*file\_path*, *recording*)

Write an annotated recording to a CSV file.

**Parameters**

- **file\_path** (*str*) – Path where the CSV file should be written.
- **recording** (*AnnotatedRecording*) – The recording data to export.

**class CountingPredicateSchmittTrigger**(*high\_predicate*, *low\_predicate*, *low\_count*=0, *high\_count*=0)

Bases: PredicateSchmittTrigger[T]

**add**(*x*)

Add new input value and return new state if changed.

**Parameters**

**x** (*T*) – Input value to evaluate.

**Returns**

New state if changed, None if no state change.

**Return type**

*SchmittState* or None

**copy\_with**(*high\_predicate*=None, *low\_predicate*=None, *low\_count*=None, *high\_count*=None)

**class CountingSchmittTrigger**(*high\_bound*, *low\_bound*, *invert*=False, *low\_count*=0, *high\_count*=0)

Bases: SchmittTrigger[T]

**add**(*x*)

Add new input value and return new state if changed.

**Parameters**

**x** (*T*) – Input value to evaluate.

**Returns**

New state if changed, None if no state change.

**Return type**

*SchmittState* or None

**copy\_with**(*high\_bound*=None, *low\_bound*=None, *low\_count*=None, *high\_count*=None, *invert*=None)

**class DelayedPredicateSchmittTrigger**(*high\_predicate*, *low\_predicate*,  
*low\_duration*=datetime.timedelta(0),  
*high\_duration*=datetime.timedelta(0))

Bases: PredicateSchmittTrigger[T]

**add(*x*)**

Add new input value and return new state if changed.

**Parameters**

***x*** (*T*) – Input value to evaluate.

**Returns**

New state if changed, None if no state change.

**Return type**

*SchmittState* or None

**copy\_with**(*high\_predicate=None, low\_predicate=None, low\_duration=None, high\_duration=None*)

**class DelayedSchmittTrigger**(*high\_bound, low\_bound, low\_duration=datetime.timedelta(0), high\_duration=datetime.timedelta(0), invert=False*)

Bases: *SchmittTrigger*[*TimedComparable*]

**add(*x*)**

Add new input value and return new state if changed.

**Parameters**

***x*** (*T*) – Input value to evaluate.

**Returns**

New state if changed, None if no state change.

**Return type**

*SchmittState* or None

**class FlexReader**

Bases: object

Universal reader that auto-detects and parses different sensor file formats.

**static parse**(*path, dist=20.0*)

Parse a sensor data file, auto-detecting the format.

Supports RSF v1, RSF v2, and CSV formats. The format is detected by examining the file header magic bytes.

**Parameters**

- **path** (*str*) – Path to the sensor data file
- **dist** (*float*) – Distance between sensor elements in mm (default: 20.0)

**Returns**

AnnotatedRecording containing the parsed measurements

**Return type**

*AnnotatedRecording*

**class MeasurementEvaluationMetric**(*value, names=None, \*, module=None, qualname=None, type=None, start=1, boundary=None*)

Bases: Enum

Enum providing standardized metrics for analyzing sensor measurements.

Each metric extracts specific features from Measurement objects, such as posture angles, movement patterns, or sensor readings.

**Metrics:**

- **LUMBAR\_ANGLE**: Measures the overall bending (flexion/extension) of the lumbar region. Calculated by summing the bottom 9 alpha values. Useful for posture and spinal movement analysis.
- **TWIST**: Measures rotational movement (twisting) of the sensor strip. Calculated by summing all beta angles. Useful for detecting torsional motions of the back.
- **LATERAL**: Measures lateral (side-to-side) flexion. Useful for analyzing sideways bending.
- **SAGITTAL**: Measures sagittal (forward/backward) flexion. Useful for analyzing bending in the sagittal plane.
- **LATERAL\_APPROX**: Provides an approximation of lateral flexion using a simplified calculation. Useful for quick or less precise assessments.
- **SAGITTAL\_APPROX**: Provides an approximation of sagittal flexion using a simplified calculation. Useful for quick or less precise assessments.
- **THORACIC\_ANGLE**: Focuses on the bending angle in the thoracic (upper back) region, typically using the second half of the sensor's angles. Useful for upper back posture analysis.
- **ACCELERATION**: Extracts the acceleration along a specific axis (usually the main axis of the sensor). Useful for detecting movement or orientation changes.
- **GYRO**: Measures angular velocity (rotation rate) from the gyroscope. Useful for analyzing dynamic movements and rotations.

**ACCELERATION = 'acceleration'**

Acceleration along main axis

**GYRO = 'gyro'**

Angular velocity (gyroscope)

**LATERAL = 'lateral'**

Lateral (side-to-side) flexion

**LATERAL\_APPROX = 'lateralApprox'**

Approximate lateral flexion

**LUMBAR\_ANGLE = 'lumbarAngle'**

Overall lumbar flexion/extension (sum of bottom 9 alpha values)

**SAGITTAL = 'sagittal'**

Sagittal (forward/backward) flexion

**SAGITTAL\_APPROX = 'sagittalApprox'**

Approximate sagittal flexion

**THORACIC\_ANGLE = 'thoracicAngle'**

Thoracic (upper back) flexion (sum of second half of alpha)

**TWIST = 'twist'**

Rotational movement (sum of all beta angles)

**exec(measurements)**

Apply this metric to a collection of measurements.

#### Parameters

**measurements** (*list* [*Measurement*] / *AnnotatedRecording*) – List of measurements or annotated recording

#### Returns

List of metric values, one per measurement

**func(*m*)**

Extract the metric value from a single measurement.

**Parameters**

**m** (*Measurement*) – The measurement to evaluate

**Returns**

The computed metric value (angle, acceleration, etc.)

**property is\_angle: bool**

**class MeasurementMetric**(*first, second*)

Bases: object

Calculates comparison metrics between two sensor measurements.

Provides various metrics to quantify differences in posture, movement, and orientation between two measurement instances.

**property alpha\_metric**

Difference in alpha angles (bending forward and backward) between measurements.

**property backward\_metric**

Difference in being bend backward

**property beta\_metric**

Difference in beta angles (rotation of the spine) between measurements.

**property forward\_metric**

Difference in being bend forward

**property lateral\_metric**

Difference between the angle of the tip compared to the case

**property metric**

Overall difference between two measurements

**property orientation\_metric**

Difference in case orientation

**property sagittal\_metric**

Difference of the angle of the sensor tip in relation to the case in the sagittal plane.

**class PredicateSchmittTrigger**(*high\_predicate, low\_predicate*)

Bases: AbstractSchmittTrigger[T]

Schmitt trigger using predicates for state transitions.

**Parameters**

- **high\_predicate** (*Callable[[T], bool]*) – Predicate to trigger transition from LOW to HIGH.
- **low\_predicate** (*Callable[[T], bool]*) – Predicate to trigger transition from HIGH to LOW.

**class RSFV1Reader**

Bases: object

Parse RSF v1 binary sensor data files into AnnotatedRecording objects.

**static parse**(*file\_path, dist=20.0*)



**class RSFV1Writer**

Bases: object

**static write**(*file\_path, recording*)

**class RSFV2Reader**

Bases: object

**static parse**(*file\_path, dist=20.0*)

**class RSFV2Writer**

Bases: object

**static write**(*file\_path, recording*)

**class SchmittState**(*value, names=None, \*, module=None, qualname=None, type=None, start=1, boundary=None*)

Bases: Enum

States for Schmitt trigger: HIGH or LOW. High indicates that filter condition is currently true.

**HIGH** = 'high'

**LOW** = 'low'

**negate**()

Return the opposite state.

**Returns**

The opposite state (HIGH <-> LOW).

**Return type**

*SchmittState*

**class SchmittTrigger**(*high\_bound, low\_bound, invert=False*)

Bases: AbstractSchmittTrigger[T]

**copy\_with**(*high\_bound=None, low\_bound=None, invert=None*)

**class TimedAnnotation**(*time, label*)

Bases: object

A time-stamped annotation label for marking events in sensor recordings.

**static from\_json\_entry**(*json*)

Create a TimedAnnotation from a JSON dictionary.

**Parameters**

**json** – Dictionary with ‘time’ and ‘value’ keys

**Returns**

TimedAnnotation instance

**Return type**

*TimedAnnotation*

**to\_json\_entry()**

Convert the annotation to a JSON-serializable dictionary.

**Returns**

Dictionary with 'time' and 'value' keys

**Return type**

dict

**class TimedComparable**(*value, time*)

Bases: Generic[T]

**create\_dataframe**(*data*)

Create a pandas DataFrame from a list of Measurement objects or an AnnotatedRecording.

**Parameters**

**data** (*List* [[Measurement](#)] or [AnnotatedRecording](#)) – The sensor data to convert. If an AnnotatedRecording, annotations will be included.

**Returns**

DataFrame with time, measurement metrics, and optional annotation column.

**Return type**

pd.DataFrame

## 9.1.1 Subpackages

### flexlib.models package

#### Subpackages

#### flexlib.models.readers package

#### Submodules

#### flexlib.models.readers.csv\_parser module

**class CSVReader**

Bases: object

Parse sensor measurement data from CSV files or strings.

**static parse**(*path, dist=20.0*)

Parse CSV content from a file.

**Parameters**

- **path** (*str*) – Path to the CSV file to parse.
- **dist** (*float*) – Distance between sensor elements in mm (default: 20.0)

**Returns**

Parsed measurements in an annotated recording.

**Return type**

[AnnotatedRecording](#)

**static parse\_string**(*content, dist=20.0*)

Parse CSV content from a string.

**Parameters**

- **content** (*str*) – The CSV content as a string.
- **dist** (*float*) – Distance between sensor elements in mm (default: 20.0)

**Returns**

Parsed measurements in an annotated recording.

**Return type**

*AnnotatedRecording*

**flexlib.models.readers.rs\_f\_v1 module****class RSFV1Reader**

Bases: object

Parse RSF v1 binary sensor data files into AnnotatedRecording objects.

**static parse**(*file\_path*, *dist*=20.0)

**flexlib.models.readers.rs\_f\_v2 module****class RSFV2Reader**

Bases: object

**static parse**(*file\_path*, *dist*=20.0)

**flexlib.models.writers package****Submodules****flexlib.models.writers.csv\_writer module****class CSVWriter**

Bases: object

A writer for exporting measurement data to CSV format.

Writes sensor values followed by timestamp in comma-separated format.

**static write**(*file\_path*, *recording*)

Write an annotated recording to a CSV file.

**Parameters**

- **file\_path** (*str*) – Path where the CSV file should be written.
- **recording** (*AnnotatedRecording*) – The recording data to export.

**flexlib.models.writers.rs\_f\_v1\_writer module****class RSFV1Writer**

Bases: object

**static write**(*file\_path*, *recording*)

**flexlib.models.writers.rs\_f\_v2\_writer module****class RSFV2Writer**

Bases: object

**static write**(*file\_path*, *recording*)

## Submodules

### flexlib.models.annotated\_recording module

**class** `AnnotatedRecording`(*measurements*, *annotations*=[])

Bases: object

A collection of sensor measurements with time-based annotations.

Contains a time-series of sensor measurements and associated annotations that mark specific events or states during the recording period.

#### **property** `end_time`

Get the timestamp of the last measurement.

**static** `from_json_entry`(*json*, *dist*=20.0)

Create an AnnotatedRecording from a JSON representation.

#### **Parameters**

- **json** – Dictionary with ‘annotations’ and ‘content’ keys
- **dist** (*float*) – Distance between sensor elements in mm (default: 20.0)

#### **Returns**

AnnotatedRecording instance

#### **Return type**

*AnnotatedRecording*

**split\_by\_annotations**()

Split the recording into segments based on annotations.

#### **Returns**

List of tuples containing (annotation\_label, measurements\_in\_segment)

#### **Return type**

*List[Tuple[str, List[Measurement]]]*

**property** `start_time`

Get the timestamp of the first measurement.

**to\_json\_entry**()

Convert the annotated recording to a JSON-serializable dictionary.

#### **Returns**

Dictionary with ‘annotations’ and ‘content’ keys

#### **Return type**

dict

**class** `TimedAnnotation`(*time*, *label*)

Bases: object

A time-stamped annotation label for marking events in sensor recordings.

**static** `from_json_entry`(*json*)

Create a TimedAnnotation from a JSON dictionary.

#### **Parameters**

- **json** – Dictionary with ‘time’ and ‘value’ keys

#### **Returns**

TimedAnnotation instance

**Return type***TimedAnnotation***to\_json\_entry()**

Convert the annotation to a JSON-serializable dictionary.

**Returns**

Dictionary with 'time' and 'value' keys

**Return type***dict***flexlib.models.flex\_reader module**

flex\_reader: Universal file reader for various sensor data formats.

Supports auto-detection and parsing of RSF v1, RSF v2, and CSV sensor files.

**class FlexReader**

Bases: object

Universal reader that auto-detects and parses different sensor file formats.

**static parse(path, dist=20.0)**

Parse a sensor data file, auto-detecting the format.

Supports RSF v1, RSF v2, and CSV formats. The format is detected by examining the file header magic bytes.

**Parameters**

- **path** (*str*) – Path to the sensor data file
- **dist** (*float*) – Distance between sensor elements in mm (default: 20.0)

**Returns**

AnnotatedRecording containing the parsed measurements

**Return type***AnnotatedRecording***flexlib.models.measurement\_evaluation\_metric module**

measurement\_evaluation\_metric: Metrics for extracting features from sensor data.

Defines the MeasurementEvaluationMetric enum for standardized feature extraction from sensor measurements.

**class MeasurementEvaluationMetric(value, names=None, \*, module=None, qualname=None, type=None, start=1, boundary=None)**

Bases: Enum

Enum providing standardized metrics for analyzing sensor measurements.

Each metric extracts specific features from Measurement objects, such as posture angles, movement patterns, or sensor readings.

**Metrics:**

- **LUMBAR\_ANGLE**: Measures the overall bending (flexion/extension) of the lumbar region. Calculated by summing the bottom 9 alpha values. Useful for posture and spinal movement analysis.
- **TWIST**: Measures rotational movement (twisting) of the sensor strip. Calculated by summing all beta angles. Useful for detecting torsional motions of the back.
- **LATERAL**: Measures lateral (side-to-side) flexion. Useful for analyzing sideways bending.
- **SAGITTAL**: Measures sagittal (forward/backward) flexion. Useful for analyzing bending in the sagittal plane.

- **LATERAL\_APPROX**: Provides an approximation of lateral flexion using a simplified calculation. Useful for quick or less precise assessments.
- **SAGITTAL\_APPROX**: Provides an approximation of sagittal flexion using a simplified calculation. Useful for quick or less precise assessments.
- **THORACIC\_ANGLE**: Focuses on the bending angle in the thoracic (upper back) region, typically using the second half of the sensor's angles. Useful for upper back posture analysis.
- **ACCELERATION**: Extracts the acceleration along a specific axis (usually the main axis of the sensor). Useful for detecting movement or orientation changes.
- **GYRO**: Measures angular velocity (rotation rate) from the gyroscope. Useful for analyzing dynamic movements and rotations.

**ACCELERATION = 'acceleration'**

Acceleration along main axis

**GYRO = 'gyro'**

Angular velocity (gyroscope)

**LATERAL = 'lateral'**

Lateral (side-to-side) flexion

**LATERAL\_APPROX = 'lateralApprox'**

Approximate lateral flexion

**LUMBAR\_ANGLE = 'lumbarAngle'**

Overall lumbar flexion/extension (sum of bottom 9 alpha values)

**SAGITTAL = 'sagittal'**

Sagittal (forward/backward) flexion

**SAGITTAL\_APPROX = 'sagittalApprox'**

Approximate sagittal flexion

**THORACIC\_ANGLE = 'thoracicAngle'**

Thoracic (upper back) flexion (sum of second half of alpha)

**TWIST = 'twist'**

Rotational movement (sum of all beta angles)

**exec(measurements)**

Apply this metric to a collection of measurements.

**Parameters**

**measurements** (*list* [[Measurement](#)] / [AnnotatedRecording](#)) – List of measurements or annotated recording

**Returns**

List of metric values, one per measurement

**func(m)**

Extract the metric value from a single measurement.

**Parameters**

**m** ([Measurement](#)) – The measurement to evaluate

**Returns**

The computed metric value (angle, acceleration, etc.)

**property is\_angle: bool**

**flexlib.models.measurement\_metric module**

measurement\_metric: Metrics for comparing two sensor measurements.

Provides the MeasurementMetric class for quantifying differences in posture, movement, and orientation. This returns a metric in the mathematical sense and thus can be used for clustering.

**class MeasurementMetric**(*first, second*)

Bases: object

Calculates comparison metrics between two sensor measurements.

Provides various metrics to quantify differences in posture, movement, and orientation between two measurement instances.

**property alpha\_metric**

Difference in alpha angles (bending forward and backward) between measurements.

**property backward\_metric**

Difference in being bend backward

**property beta\_metric**

Difference in beta angles (rotation of the spine) between measurements.

**property forward\_metric**

Difference in being bend forward

**property lateral\_metric**

Difference between the angle of the tip compared to the case

**property metric**

Overall difference between two measurements

**property orientation\_metric**

Difference in case orientation

**property sagittal\_metric**

Difference of the angle of the sensor tip in relation to the case in the sagittal plane.

**flexlib.models.schmitt\_trigger module**

schmitt\_trigger: Classes for robust edge detection, debouncing, and filtering using Schmitt triggers.

Includes: - AbstractSchmittTrigger: Base class for Schmitt triggers - PredicateSchmittTrigger: Predicate-based triggers - CountingPredicateSchmittTrigger: Adds count-based noise immunity - DelayedPredicateSchmittTrigger: Adds time-based filtering - SchmittTrigger: Classic numeric Schmitt trigger - CountingSchmittTrigger: Numeric + count-based filter - DelayedSchmittTrigger: Numeric + time-based filter

These classes are useful for signal processing, user input handling, and control systems.

**class AbstractSchmittTrigger**

Bases: Generic[T]

Base class for Schmitt triggers with hysteresis-based state detection. The Schmitt trigger classes implementations for robust signal state detection based on hysteresis. Schmitt triggers are fundamental components in signal processing and digital electronics, used to convert noisy or analog input signals into clean, discrete digital outputs. They achieve this by introducing hysteresis, which means the trigger has two different threshold levels: one for transitioning from LOW to HIGH (the upper threshold), and another for transitioning from HIGH to LOW (the lower threshold). This prevents rapid toggling or “chattering” of the output when the input signal is noisy or hovers near the threshold. This module provides several flexible and extensible Schmitt trigger implementations:

**\_state**

The current state of the trigger (HIGH or LOW).

**Type**

SchmittState

**add(*x*)**

Add new input value and return new state if changed.

**Parameters**

**x** (*T*) – Input value to evaluate.

**Returns**

New state if changed, None if no state change.

**Return type**

SchmittState or None

**property state: SchmittState**

Current trigger state.

**Returns**

The current state (HIGH or LOW).

**Return type**

SchmittState

**class CountingPredicateSchmittTrigger**(*high\_predicate*, *low\_predicate*, *low\_count*=0, *high\_count*=0)

Bases: PredicateSchmittTrigger[T]

**add(*x*)**

Add new input value and return new state if changed.

**Parameters**

**x** (*T*) – Input value to evaluate.

**Returns**

New state if changed, None if no state change.

**Return type**

SchmittState or None

**copy\_with**(*high\_predicate*=None, *low\_predicate*=None, *low\_count*=None, *high\_count*=None)

**class CountingSchmittTrigger**(*high\_bound*, *low\_bound*, *invert*=False, *low\_count*=0, *high\_count*=0)

Bases: SchmittTrigger[T]

**add(*x*)**

Add new input value and return new state if changed.

**Parameters**

**x** (*T*) – Input value to evaluate.

**Returns**

New state if changed, None if no state change.

**Return type**

SchmittState or None

**copy\_with**(*high\_bound*=None, *low\_bound*=None, *low\_count*=None, *high\_count*=None, *invert*=None)



```
class DelayedPredicateSchmittTrigger(high_predicate, low_predicate,  
                                     low_duration=datetime.timedelta(0),  
                                     high_duration=datetime.timedelta(0))
```

Bases: PredicateSchmittTrigger[T]

**add**(*x*)

Add new input value and return new state if changed.

**Parameters**

**x** (T) – Input value to evaluate.

**Returns**

New state if changed, None if no state change.

**Return type**

SchmittState or None

**copy\_with**(*high\_predicate=None, low\_predicate=None, low\_duration=None, high\_duration=None*)

```
class DelayedSchmittTrigger(high_bound, low_bound, low_duration=datetime.timedelta(0),  
                             high_duration=datetime.timedelta(0), invert=False)
```

Bases: SchmittTrigger[TimedComparable]

**add**(*x*)

Add new input value and return new state if changed.

**Parameters**

**x** (T) – Input value to evaluate.

**Returns**

New state if changed, None if no state change.

**Return type**

SchmittState or None

```
class HasDateTime(*args, **kwargs)
```

Bases: Protocol

**datetime**: datetime

```
class PredicateSchmittTrigger(high_predicate, low_predicate)
```

Bases: AbstractSchmittTrigger[T]

Schmitt trigger using predicates for state transitions.

**Parameters**

- **high\_predicate** (Callable[[T], bool]) – Predicate to trigger transition from LOW to HIGH.
- **low\_predicate** (Callable[[T], bool]) – Predicate to trigger transition from HIGH to LOW.

```
class SchmittState(value, names=None, *, module=None, qualname=None, type=None, start=1,  
                  boundary=None)
```

Bases: Enum

States for Schmitt trigger: HIGH or LOW. High indicates that filter condition is currently true.

**HIGH** = 'high'

**LOW** = 'low'

**negate()**

Return the opposite state.

**Returns**

The opposite state (HIGH <-> LOW).

**Return type**

SchmittState

**class SchmittTrigger**(*high\_bound, low\_bound, invert=False*)

Bases: AbstractSchmittTrigger[T]

**copy\_with**(*high\_bound=None, low\_bound=None, invert=None*)

**class TimedComparable**(*value, time*)

Bases: Generic[T]

**get\_schmitt\_episodes**(*trigger*)

## 9.1.2 Submodules

### 9.1.3 flexlib.dataframe\_builder module

dataframe\_builder: Utilities for converting sensor data to pandas DataFrames for analysis.

**create\_dataframe**(*data*)

Create a pandas DataFrame from a list of Measurement objects or an AnnotatedRecording.

**Parameters**

**data** (*List* [Measurement] or AnnotatedRecording) – The sensor data to convert. If an AnnotatedRecording, annotations will be included.

**Returns**

DataFrame with time, measurement metrics, and optional annotation column.

**Return type**

pd.DataFrame

## DEMO: PLOTTING

This notebook contains a few simple plots that can be created with flexlib

```
[2]: from flexlib import Measurement, FlexReader
dist = 25
data = FlexReader().parse("../docs/test_data/tidy_up.rs", dist = dist)
m = data[0]
```

This create a simple 2D of the sensor viewed from the right.

```
[3]: import matplotlib.pyplot as plt
from mpl_toolkits.mplot3d import Axes3D


max_val = 18 * dist

# Extract coordinates
# Plot 3D coordinates
fig = plt.figure(figsize=(8, 6))
ax = fig.add_subplot(111)

ax.plot(m.coordinates.y, m.coordinates.z, marker='o')

for line in ax.get_lines():
    line.set_linewidth(3)
ax.set_xlabel('Y (mm)')
ax.set_ylabel('Z (mm)')
ax.set_title('3D Sensor Coordinates')
ax.set_xlim(-max_val, max_val)
ax.set_ylim(0, max_val)
```

```
[3]: (0.0, 450.0)
```



```
[4]: def equally_spaced_sample(lst, n):
      if n >= len(lst):
          return lst
      idxs = [round(i * (len(lst) - 1) / (n - 1)) for i in range(n)]
      return [lst[i] for i in idxs]
```

This as a kind of 2D-Heatmap with 500 sensor positions overlayed on each other.

```
[5]: import matplotlib.pyplot as plt
      from mpl_toolkits.mplot3d import Axes3D

      max_val = 18 * dist

      # Extract coordinates
      # Plot 3D coordinates
      fig = plt.figure(figsize=(8, 6))
      ax = fig.add_subplot(111)

      for m in equally_spaced_sample(data.measurements, 500):
          ax.plot(m.coordinates.y, m.coordinates.z, color="#00000001")


      for line in ax.get_lines():
          line.set_linewidth(10)
      ax.set_xlabel('Y (mm)')
      ax.set_ylabel('Z (mm)')
      ax.set_title('3D Sensor Coordinates')
```

(continues on next page)

(continued from previous page)

```
ax.set_xlim(-max_val, max_val)
ax.set_ylim(0, max_val)
```

```
[5]: (0.0, 450.0)
```




This is a heatmap of the person viewed from above. In this case the person was primarily bend forward and head more movement in the sagittal plane than in the lateral.

```
[6]: import numpy as np
import seaborn as sns

from math import pi

f = 180 / pi

# Concatenate all y and z coordinates from all measurements
all_y = [f * m.lateral_flexion for m in data.measurements]
all_z = [f * m.sagittal_flexion for m in data.measurements]
sns.jointplot(x=all_y, y=all_z, kind="hist", marginal_kws=dict(bins=30, fill=True))
plt.xlabel('Lateral (Y, mm)')
plt.ylabel('Sagittal (Z, mm)')
plt.suptitle('Marginal Histogram of Lateral and Sagittal Values', y=1.02)
plt.xlim(-120, 120)
plt.ylim(-120, 120)
plt.show()
```



```
[7]: def measurement_to_3d_coords(m : Measurement):
    reconstruction = m.full_reconstruction

    # Stack left and right edge coordinates into arrays
    left = np.vstack([reconstruction.left.x, reconstruction.left.y, reconstruction.
↪ left.z]).T
    right = np.vstack([reconstruction.right.x, reconstruction.right.y, reconstruction.
↪ right.z]).T

    return left, right, np.vstack([left[:, 0], right[:, 0]]), np.vstack([left[:, 1],
↪ right[:, 1]]), np.vstack([left[:, 2], right[:, 2]])
```

This code creates a 3D-Plot of the sensor

```
[8]: import matplotlib.pyplot as plt
    from mpl_toolkits.mplot3d import Axes3D
    import numpy as np
```

(continues on next page)

(continued from previous page)

```
# Use the full reconstruction object


fig = plt.figure(figsize=(10, 7))
ax = fig.add_subplot(111, projection='3d')

left, right, X, Y, Z = measurement_to_3d_coords(m)
# Plot the 3D foil as a surface
ax.plot_surface(X, Y, Z, color='gray', alpha=0.5, rstride=1, cstride=1, linewidth=0)

# Optionally plot the left and right edges
ax.plot(left[:, 0], left[:, 1], left[:, 2], color='gray', alpha=1, linewidth=2)
ax.plot(right[:, 0], right[:, 1], right[:, 2], color='gray', alpha=1, linewidth=2)

ax.set_xlabel('X (mm)')
ax.set_ylabel('Y (mm)')
ax.set_zlabel('Z (mm)')

lim = dist * 18
ax.set_xlim(-lim, lim)
ax.set_ylim(-lim, lim)
ax.set_zlim(0, lim)
ax.set_title('3D Flextail Full Reconstruction')
plt.tight_layout()
plt.show()
```





## DISHWASHER DEMO

This notebook provides a straightforward overview of ergonomic analysis for unloading a dishwasher, using sensor data to detect bending periods and evaluate posture quality. The data features here is a person that is unloading a dishwasher for 1 minute measured with a frequency of 25 Hz.

The data analysis includes:

- Loading and processing time series data for body angles and movement.
- Detecting “bending forward” actions
- Scoring posture quality during each bending period based on lumbar angle and torso twist.

```
[1]: # flexlib is a python library written by Minktec to facilitate easy flexlib data processing
import flexlib
from flexlib import MeasurementEvaluationMetric as mem
from flexlib import CountingPredicateSchmittTrigger, SchmittState

import seaborn as sns # plotting
import matplotlib.pyplot as plt # plotting
import pandas as pd # dataframes / tables
import numpy as np # math
```

After importing the libraries, we can load the data recorded with the FlexTail-Sensor. The `create_dataframe` method takes the raw values and computes metrics derived from the sensor data. This makes it easy to work with the sensor data. If more complicated evaluations are necessary, the sensor angles or coordinates can also be computed.

```
[2]: measurements : flexlib.AnnotatedRecording = flexlib.FlexReader().parse("../docs/test_
data/dishwasher.rsfs")
df : pd.DataFrame = flexlib.create_dataframe(measurements)
df
```

```
e:\flexlib\.venv\Lib\site-packages\flexlib\dataframe_builder.py:47: FutureWarning:
Setting an item of incompatible dtype is deprecated and will raise an error in a
future version of pandas. Value 'loading dishwasher' has dtype incompatible with
float64, please explicitly cast to a compatible dtype first.
df.loc[closest_idx, 'annotation'] = annotation.label
```

```
[2]:
```

	Time	annotation	lumbarAngle	twist	\
0	2024-12-10 08:30:19.051	loading dishwasher	-0.302237	-0.032632	
1	2024-12-10 08:30:19.118	NaN	-0.306199	-0.017651	
2	2024-12-10 08:30:19.169	NaN	-0.305224	-0.021672	
3	2024-12-10 08:30:19.200	NaN	-0.279675	0.032262	
4	2024-12-10 08:30:19.258	NaN	-0.318241	-0.013603	
...	...	...	...	...	
1480	2024-12-10 08:31:19.709	NaN	0.119518	0.103478	
1481	2024-12-10 08:31:19.743	NaN	0.123758	0.083975	
1482	2024-12-10 08:31:19.829	NaN	0.125396	0.081779	

(continues on next page)

(continued from previous page)

```

1483 2024-12-10 08:31:19.860      NaN      0.124641  0.091700
1484 2024-12-10 08:31:19.888      NaN      0.123758  0.086738

      lateral  sagittal  lateralApprox  sagittalApprox  thoracicAngle  \
0      0.071393 -0.064038      0.046761      -0.027765      0.310219
1      0.062159 -0.061405      0.034435      -0.025203      0.306267
2      0.070013 -0.053050      0.041685      -0.016517      0.308264
3      0.101480 -0.024342      0.069655      0.007804      0.311819
4      0.077045 -0.064217      0.044471      -0.031815      0.311717
...      ...      ...      ...      ...      ...
1480 0.370382  0.748794      0.340798      0.833886      0.474724
1481 0.373323  0.773125      0.341342      0.858399      0.474082
1482 0.402006  0.767567      0.369646      0.852625      0.469927
1483 0.416074  0.769438      0.381614      0.854627      0.468836
1484 0.407534  0.777981      0.374646      0.863071      0.471117

      acceleration  gyro
0      0.9830  1.66
1      1.0015  1.36
2      0.9830  1.16
3      0.9675  0.80
4      1.0000  0.26
...      ...      ...
1480 0.8685  3.14
1481 0.8670  3.46
1482 0.8530  2.94
1483 0.8560  2.84
1484 0.8405  2.06

[1485 rows x 11 columns]

```

**Lumbar Angle:**

The lumbar angle is computed from sensor data placed on the lower back, representing the flexion or extension of the lumbar. It quantifies how much the lower back bends forward or backward during movement.

**Sagittal**


The sagittal angle is the value of the torso's forward/backward tilt in the sagittal plane, derived from the sensor shape and the case position.

Both values are reported in radians where positive means bending forward.

```

[3]: plt.figure(figsize=(12, 6))
plt.plot(df['Time'], df['lumbarAngle'], label='Lumbar Angle')
plt.plot(df['Time'], df['sagittal'], label='Sagittal')
plt.xlabel('Time')
plt.ylabel('Angle (radians)')
plt.title('Lumbar Angle and Sagittal Over Time')
plt.legend()
plt.tight_layout()
plt.show()

```



## 11.1 Schmitt Trigger for Bending Detection

Next, we want to know when the person is bending forward. For this, we use a [Schmitt trigger](#).

This allows us to efficiently filter out noise or jitter caused by the person.

If we were to count everything as bending where the sagittal angle is  $> 30$  degrees, it could lead to very short or undefined edges in the detected timestamps. The Schmitt trigger allows us to set a trigger threshold (high), in this example 30 degrees, and a release threshold (low), for example 25 degrees.

In this case, we use a `CountingSchmittTrigger`, which only changes the state when the set thresholds are met for a certain number of measurements.

[7]:

```
# Define thresholds and counts for the Schmitt trigger
sagittal_high = 0.15
sagittal_low = 0.05
high_count = 3 # Number of consecutive points above high to trigger
low_count = 3 # Number of consecutive points below low to reset

# Define predicates for the sagittal signal
high_predicate = lambda x: x > sagittal_high
low_predicate = lambda x: x < sagittal_low

# Create the CountingPredicateSchmittTrigger
trigger = CountingPredicateSchmittTrigger(
    high_predicate=high_predicate,
    low_predicate=low_predicate,
    high_count=high_count,
    low_count=low_count
)

# Track the trigger state for each time point
bending_mask = []
for val in df['sagittal']:
    trigger.add(val)
```

(continues on next page)

(continued from previous page)


```

bending_mask.append(trigger.state == SchmittState.HIGH)

df['bending'] = bending_mask

plt.figure(figsize=(12, 6))
plt.plot(df['Time'], df['lumbarAngle'], label='Lumbar Angle')
plt.plot(df['Time'], df['sagittal'], label='Sagittal')
plt.fill_between(df['Time'], df['lumbarAngle'].min(), df['sagittal'].max(), where=df[
    ↪ 'bending'], color='orange', alpha=0.2, label='Bending Forward')
plt.xlabel('Time')
plt.ylabel('Sagittal Angle (radians)')
plt.title('Sagittal Angle with Bending Forward Periods')
plt.legend(loc='upper left')
plt.tight_layout()
plt.show()

```



Great. We can clearly see that the periods marked as being bend forward match our expectations and are clearly separated. In this case one “Bending Forward” block represents taking a mug or plate out of the dishwasher.

Next we want to do a simple evaluation of the ergonomic quality of the movements. For that we use a combination of the bending of the lumbar spine and the rotation of the thorso twist. We are going to rate the movements from 0 (bad) to 1 (good).

```

[9]: lumbar = df['lumbarAngle']
     twist = df['twist']

     # Score: 1 is perfect, 0 is bad; penalize more when both are bad
     # This score is kind of arbitrary and can be adjusted based on the specific_
     ↪ requirements
     lumbar_norm = (lumbar / 0.5).clip(0, 1)
     twist_norm = (twist.abs() / 0.15).clip(0, 1)
     df['score'] = 1 - (lumbar_norm + twist_norm) / 3 - (lumbar_norm * twist_norm)

     # Find contiguous bending forward periods using pandas magic
     bending_periods = [group.index.tolist() for _, group in df[df['bending']].groupby((~
     ↪ df['bending']).cumsum())]

```

(continues on next page)

(continued from previous page)


```
# Compute average score for each bending period (simple version)
period_scores = []
for i, period in enumerate(bending_periods):
    avg = df['score'].loc[period].mean()
    period_scores.append({'period': str(i+1), 'score': avg})

period_scores_df = pd.DataFrame(period_scores)

plt.figure(figsize=(10, 5))
plt.bar(period_scores_df['period'], period_scores_df['score'])
plt.xlabel('Bending Forward Period')
plt.ylabel('Average Score')
plt.title('Bending Forward Quality per Period')
plt.ylim(-0.05, 1.05)

x_pos = 5.5
plt.axvline(x=x_pos)
plt.text(x_pos+0.05, 0.95, 'Start emptying bottom drawer', rotation=90, va='top',
        ↪ fontsize=11, fontweight='bold')

plt.tight_layout()
plt.show()
```



Based on our simple score we can clearly see that once the person starts unloading the bottom tray of the dishwasher, the ergonomic quality starts to degrade.

## 11.2 Summary

We have:

1. Loaded sensor data
2. Created a plot
3. Detected bending forward periods
4. Evaluated the ergonomic impact